

Beer's Law Calibration Curve

Because crystal violet is a colored species, a device that measures light absorbance may be used to determine concentration. The absorbance is related to concentration according to the Beer-Lambert law:

$$A = \epsilon lc$$

where A is the absorbance, ϵ is the “molar absorptivity” coefficient (a constant unique to the solute at a particular wavelength) for the substance of interest, l is the length of the light path, and c is the concentration of the absorbing species. A plot of absorbance vs. concentration for a set of solutions will yield a straight line for which the slope is ϵl . The cells are 1 cm wide, so the l -term is ignored. This plot can then be used to find concentration given absorbencies.

To prepare a Beer's law calibration curve, you should first determine the wavelength of light to use for maximum absorption. Then determine the absorbencies for a set of at least seven different concentrations of crystal violet prepared by serially diluting the stock solution – 5.00 mL CV⁺ solution with 5.00 mL water for a total volume of 10.00 mL.

Graph your results to determine the linear range between concentration and absorbance.

Wavelength: _____

Molarity stock: _____

| Solution Number | Concentration (M) of Crystal Violet | Absorbance |
|-----------------|---|------------|
| 1 | | |
| 2 | | |
| 3 | | |
| 4 | | |
| 5 | | |
| 6 | | |
| 7 | | |
| 8 | | |

Linear regression equation for absorbance vs. concentration:

Correlation coefficient and reason for discarding data: